

DLLComm

Dynamic linked library




What is DLLComm?

Instrument for connecting of user's datalogging software with control system of ovens and incubators (except CO2Cell) allowing to continuously control and receive data from the MMM chamber.

How does it work?

- 1) Purchase the DLLComm software packet (order code MC 001980)
- 2) Install the software according to user's software environment
- 3) Define data connection between a user's data logging software and MMM device in source code (C++ skills needed)



Main features

- Easy instrument for controlling and receiving data within C++ environment
- Possibility of creating an interface between your information system and connected devices from our production
- Possibility to control the device, including programs settings, starting, stopping, recording, etc.
- Recorded data can be easily integrated into your IT systems and it is up to your consideration, how else to use the data
- You can manage the users policies and other advanced settings, too
- Comprehensive user manual

DLLComm is currently available for Windows only.
If you want to use it on different operating system, please write a request to us.



www.mmm-medcenter.de

MMM Medcenter Einrichtungen GmbH
Semmelweisstrasse 6 D-82152 Planegg / München
Tel.: +49 89 8992 2620
E-mail: medcenter@mmm-medcenter.de
www.mmm-medcenter.de